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MEDICAL COMPANIES

The Function of Surgical Stapling

- Primary Objective: Hemostasis and leak prevention
- Cause of Leak and Hemorrhage
 - Wrong staple height selection (cartridge or reload)

Tissue Dynamics

Two physical states of living tissue must be considered: the fluid phase (extra cellular water is present in the tissue prior to compression); and the solid phase (after fluid evacuation, the tissue behaves much like a solid)

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HYSTERESIS

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There is a time lag between when pressure is applied to tissue and when the tissue's elastic forces can achieve balance with this pressure

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
Johnson & Johnson MEDICAL DEVICES COMPANIES

What Tissue like?


Johnson & Johnson MEDICAL DEVICES COMPANIES

EES Approach-Compression


Unmet Need in Stapling: Tissue Slippage

Staple Formation

Tissue Slippage


GST System design enhancements

7

Gripping Surface Technology

Improved staple trajectory to the pocket by:

- ▶ Reducing the impact of tissue flow
- ▶ Support staple legs for a longer distance
- ▶ Improved longitudinal and lateral alignment

Asymmetric, bent tips

improve the tip trajectory to the pockets


ETHICON

PART OF THE Johnson & Johnson FAMILY OF COMPANIES

8

Tissue slippage effect on staple formation

Providing alignment for staple legs to hit the anvil pockets


Tissue flow demonstration during firing with gold reload at tissue thickness of 3.2mm

9

Less Tissue Slippage for More Targeted Transections

Less tissue milking to potentially eliminate reloads[‡]

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
johnson-johnson MEDICAL DEVICES COMPANIES

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Evaluation of a Powered Stapler System with Gripping Surface Technology on Surgical Interventions Required During Laparoscopic Sleeve Gastrectomy.

Fegelman E¹, Knippenberg S¹, Schwiers M¹, Stefanidis D², Gersin KS², Scott JD³, Fernandez AZ⁴.

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
Preclinical quantification of air leaks in a physiologic lung model: effects of ventilation modality and staple design.

Eckert CE¹, Harris JL¹, Wong JB¹, Thompson S², Kassis ES³, Tsuboi M⁴, Ott HC⁵, Force S⁶.

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Notes: ^aMann-Whitney test. ^bTwo-sample proportion test. ^cANOVA test. ^dt-test. Mean (SD) leak volume, leak rate, leak incidence rate, and trachea pressure for PPV testing. P-values from statistical tests are provided. Note: in eight of the datasets (four GRD, four UNI specimens), transducer calibration errors resulted in corrupt trachea pressures, so they were not included in the trachea pressure statistical testing.

Abbreviations: GRD, graduated staple height; PPV, positive pressure ventilation; UNI, uniform staple height.

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Table 3 Summary leak data for NPV

Test group	Leak volume (mL)	Leak rate (mL/min)	P-value	Leak incidence rate (%)	P-value	Chamber pressure (cm H ₂ O)	P-value
GRD (thinner)	57.1 (202)	5.71 (20.2)	0.014 ^a	44	0.015 ^b	-14.1 (3.63)	0.817 ^c
UNI (thinner)	2.02 (5.57)	0.202 (0.557)		16		-13.9 (2.79)	
GRD (thicker)	242 (711)	24.2 (71.1)		43		-14.7 (2.54)	
UNI (thicker)	32.7 (86.6)	3.27 (8.66)		23		-14.5 (2.50)	
GRD (thinner + thicker)	158 (546)	15.8 (54.6)	0.003 ^a	44	0.004 ^b	-14.3 (3.07)	0.858 ^d
UNI (thinner + thicker)	18.7 (65.4)	1.87 (6.54)		20		-14.3 (2.62)	

Notes: ^aMann-Whitney test. ^btwo-sample proportion test. ^cANOVA test. ^dt-test. Mean (SD) leak volume, leak rate, leak incidence rate, and chamber pressure for NPV testing. P-values from statistical tests are provided.

Abbreviations: GRD, graduated staple height; NPV, negative pressure ventilation; UNI, uniform staple height.

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Johnson & Johnson MEDICAL DEVICES COMPANIES

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Original Article - Research
Veterinary Surgery / Volume 45, Issue 1 First published: 23 November 2015

Inflation Pressures for Ex Vivo Lung Biopsies After Application of Graduated Compression Staples

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Obes Surg. 2018 Nov;28(11):3446-3453. doi: 10.1007/s11695-018-3363-0.

Initial Assessment of Mucosal Capture and Leak Pressure After Gastrointestinal Stapling in a Porcine Model.

Thompson SE¹, Young MT¹, Lewis MT¹, Boronyak SM¹, Clymer JW², Fegelman EJ¹, Nagle DA¹.

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